



Resin-based root filling material

An easy, accurate root filling system for an **improved long-lasting seal**

Description

BioRaCe Fill comprises a resin-based core material that will bond to the resin sealer that in turn bonds to clean dentin leaving virtually no room for leakage.

The new BioRaCe Fill Self-Etch Sealer* eliminates the priming step without hindering the sealer's performance.

And when switching from gutta-percha to the BioRaCe Fill, there's absolutely no change in technique needed.

BioRaCe Fill offers the following features:

- Biocompatibility: nontoxic and nonmutagenic (see available online literature)
- Improved leak-resistance: unlike gutta-percha, BioRaCe Fill leaves no consistent gap for leakage
- Technique-Compatible: works with your current filling method
- Superior radiopacity
- Retreatmentability: easy with solvent and e.g. D-RaCe desobturation instruments

Availability

- Syringes of Self-Etch Sealer (4ml)
- Standardised Points #25/.06 - #35/.02 - #35/.04 - #40/.02 - #40/.04 - #50/.04 - #60/.02 (50 pces/box)
- Accessory points FM (50 pces/box)
- Mixing tips

An Intro Kit is available! It contains:

- 1 syringe of Self-Etch Sealer (4ml);
- 1 box of standardised points #35/.04 and 1 of #40/.04 (50 pces/box);
- 1 box of standardised assorted points (80 pces/box):
#35/.02 - #40/.02 - #50/.04 - #60/.02;
- 1 box of accessory points FM (120 pces/box);
- 20 mixing tips;
- 1 instruction manual.



Website

www.biorace.ch

* Primary component of BioRaCe Fill root filling system: Resilon II™